


Anti-Vegf Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P15692

Item Number: CM0004342



Introduction

Anti-VEGF rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat VEGF (15-27/34-45 kDa). Ideal for angiogenesis and hypoxia vascular research.

[Learn More](#)

Parameter	Value
Product Name	VEGF Rabbit mAb
Remarks/Alias	VEGFA, MVCD1, VEGF, VPF, vascular endothelial growth factor A
Species	Human
GeneID (Human)	7422
GeneID	7422
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-135 of human VEGF (NP_001165094.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15692
Protein Weight	15-27kDa/34-45kDa
Shipping	Ice bag