

Anti-Vegf Receptor 2 Polyclonal Antibody For Wb, Elisa - P35968

Item Number: CM0014829



Introduction

Rabbit polyclonal antibody to VEGF Receptor 2 (KDR/FLK1/CD309) validated for Western blot and ELISA. Detects human, mouse, rat. Immunogen: recombinant fusion protein (aa 1147-1356). Ideal for angiogenesis and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	VEGF Receptor 2 Rabbit pAb
Remarks/Alias	FLK1; CD309; VEGFR; VEGFR2; VEGF Receptor 2
Species	Human
GeneID (Human)	3791
GeneID	3791
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1147-1356 of human VEGF Receptor 2 (NP_002244.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35968
Protein Weight	152kDa
Shipping	Ice bag