


Anti-Vegfa Mouse Monoclonal Antibody For Wb, Elisa, Ihc-P - P15692

Item Number: CM0002783



Introduction

High-quality mouse monoclonal antibody against human VEGFA, validated for Western blot, ELISA, and immunohistochemistry with cross-reactivity to mouse and rat. Ideal for angiogenesis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	VEGFA Mouse mAb
Remarks/Alias	VPF; VEGF; MVCD1; VEGFA
Species	Human
Gene ID	7422
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 207-327 of human VEGFA (P15692).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA, IHC-P
Cross Reactivity	Mouse, Rat
SWISS	P15692
Protein Weight	27kDa
Shipping	Ice bag