


Anti-Vegfa Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15692

Item Number: CM0018406



Introduction

Rabbit polyclonal antibody against human VEGFA, suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for angiogenesis and hypoxia research.

[Learn More](#)

Parameter	Value
Product Name	VEGFA Rabbit pAb
Remarks/Alias	VPF; VEGF; MVCD1; VEGFA
Species	Human
GeneID (Human)	7422
GeneID	7422
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human VEGFA (NP_001165094.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15692
Protein Weight	27kDa
Shipping	Ice bag