

Anti-Vegfb Polyclonal Antibody For Wb, If/Icc, Elisa - P49765

Item Number: CM0020978



Introduction

Rabbit polyclonal antibody against human VEGFB. Suitable for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Targets vascular endothelial growth factor B, a key regulator of angiogenesis and endothelial cell function.

[Learn More](#)

Parameter	Value
Product Name	VEGFB Rabbit pAb
Remarks/Alias	VRF; VEGFL; VEGFB
Species	Human
Gene ID (Human)	7423
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-121 of human VEGFB (NP_003368.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49765
Protein Weight	22kDa
Shipping	Ice bag