

Anti-Vegfb Polyclonal Antibody For Wb, Ihc-P, Elisa - P49765

Item Number: CM0015955



Introduction

Rabbit polyclonal antibody against human VEGFB (VRF, VEGFL). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for angiogenesis and cardiovascular studies.

[Learn More](#)

Parameter	Value
Product Name	VEGFB Rabbit pAb
Remarks/Alias	VRF; VEGFL; VEGFB
Species	Human
Gene ID (Human)	7423
Gene ID	7423
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse VEGFB (NP_035827.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49765
Protein Weight	22kDa
Shipping	Ice bag