


Anti-Vegfc (Flt4-L, Vrp) Rabbit Monoclonal Antibody For Wb And Elisa - P49767

Item Number: CM0006504



Introduction

CamelBio VEGFC Rabbit mAb is a high-specificity monoclonal antibody for WB and ELISA. Recognizes human VEGFC within amino acids 112-227. Ideal for angiogenesis, lymphangiogenesis, and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	VEGFC Rabbit mAb
Remarks/Alias	VRP; Flt4-L; LMPH1D; LMPHM4; VEGFC
Species	Human
GeneID (Human)	7424
GeneID	7424
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 112-227 of human VEGFC. (NP_005420.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49767
Protein Weight	47kDa
Shipping	Ice bag