


Anti-Vegfc Rabbit Polyclonal Antibody For Wb, Elisa - P49767

Item Number: CM0025778



Introduction

Rabbit polyclonal anti-VEGFC antibody for Western blot and ELISA. Detects human, mouse, rat VEGFC (47 kDa). Ideal for angiogenesis and lymphangiogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	VEGFC Rabbit pAb
Remarks/Alias	VRP; Flt4-L; LMPH1D; LMPHM4; VEGFC
Species	Human
Gene ID (Human)	7424
Gene ID	7424
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human VEGFC (NP_005420.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49767
Protein Weight	47kDa
Shipping	Ice bag