

Anti-Vegfr1 Rabbit Monoclonal Antibody For Wb, Elisa - P17948

Item Number: CM0002709



Introduction

Rabbit monoclonal antibody against human VEGFR1 (FLT-1), validated for WB and ELISA. Cross-reacts with mouse and rat. For angiogenesis, cancer, and developmental research.

[Learn More](#)

Parameter	Value
Product Name	VEGFR1 Rabbit mAb
Remarks/Alias	FLT; FLT-1; VEGFR1; VEGFR-1; FRT
Species	Human
GeneID (Human)	2321
GeneID	2321
Immunogen	Synthetic Peptide. A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VEGFR1 (P17948).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P17948
Protein Weight	151kDa
Shipping	Ice bag