


Anti-Vegfr1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17948

Item Number: CM0013031



Introduction

High-quality rabbit polyclonal antibody against human VEGFR1 (FLT-1), suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for angiogenesis, cancer, and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	VEGFR1 Rabbit pAb
Remarks/Alias	FLT; FLT-1; VEGFR1; VEGFR-1
Species	Human
GeneID (Human)	2321
GeneID	2321
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human VEGFR1 (NP_002010.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17948
Protein Weight	151kDa
Shipping	Ice bag