

Anti-Vegfr2/Cd309 Monoclonal Antibody For Fc - P35968

Item Number: CM0008114



Introduction

Fluorescent () rabbit monoclonal antibody against human VEGFR2/CD309 (KDR/FLK1), validated for flow cytometry (FC). Detects receptor tyrosine kinase involved in angiogenesis and vascular development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human VEGFR2/CD309 mAb
Remarks/Alias	FLK1; CD309; VEGFR; VEGFR2
Species	Human
Gene ID (Human)	3791
Gene ID	3791
Immunogen	Recombinant protein of human VEGFR2/CD309
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P35968
Protein Weight	152kDa
Shipping	Ice bag