

Anti-Vegfr2/Kdr/Flk-1/Cd309 Rabbit Polyclonal Antibody For Wb And Elisa - P35918

Item Number: CM0022151

Introduction

Anti-VEGFR2/KDR/Flk-1/CD309 Rabbit Polyclonal Antibody. Validated for Western Blot and ELISA with mouse samples. Targets vascular endothelial growth factor receptor 2, crucial for angiogenesis research. High-quality reagent for studying angiogenesis, vascular development, and VEGFR signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	VEGFR2/KDR/Flk-1/CD309 Rabbit pAb
Remarks/Alias	orv; Flk1; Ly73; Flk-1; Krd-1; VEGFR2; VEGFR-2; sVEGFR-2; 6130401C07; VEGFR2/KDR/Flk-1/CD309
Species	Mouse
Gene ID (Human)	3791
Gene ID	16542
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-762 of mouse VEGFR2/KDR/Flk-1/CD309 (NP_034742.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P35918
Protein Weight	150kDa
Shipping	Ice bag