


Anti-Vegfr3/Flt4 Rabbit Polyclonal Antibody For Wb And Elisa - P35916

Item Number: CM0016123



Introduction

Rabbit polyclonal antibody against human VEGFR3/FLT4, validated for WB and ELISA. Targets a synthetic peptide within amino acids 1230-1329. Suitable for lymphangiogenesis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	VEGFR3/FLT4 Rabbit pAb
Remarks/Alias	PCL; CHTD7; FLT-4; FLT41; LMPH1A; LMPHM1; VEGFR3; VEGFR-3; VEGFR3/FLT4
Species	Human
Gene ID (Human)	2324
Gene ID	2324
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1230-1329 of human VEGFR3/FLT4 (NP_891555.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35916
Protein Weight	153kDa
Shipping	Ice bag