

Anti-Ventx Polyclonal Antibody For Wb, Elisa - O95231

Item Number: CM0026937

Introduction

Rabbit polyclonal antibody against human VENTX, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for developmental biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	VENTX Rabbit pAb
Remarks/Alias	NA88A, HPX42B, VENTX2, VENTX
Species	Human
Gene ID (Human)	27287
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-258 of human VENTX (NP_055283.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95231
Protein Weight	28kDa
Shipping	Ice bag