

Anti-Veph1 Rabbit Polyclonal Antibody For Wb, Elisa - Q14D04

Item Number: CM0031080

Introduction

VEPH1 Rabbit Polyclonal Antibody validated for WB, ELISA. Reacts with human and mouse VEPH1, a membrane protein inhibiting TGF-beta, FOXO, Hippo, and Wnt signaling. Ideal for IVD research and bulk supply.

[Learn More](#)

Parameter	Value
Product Name	VEPH1 Rabbit pAb
Remarks/Alias	MELT; VEPH; VEPH1
Species	Human
GeneID (human)	79674
GeneID	79674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 448-547 of human VEPH1 (NP_078897.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14D04
Protein Weight	95kDa
Shipping	Ice bag