

Anti-Versican Monoclonal Antibody For Wb, Elisa - P13611

Item Number: CM0002539

Introduction

Rabbit monoclonal antibody against human Versican (CSPG2). Suitable for WB and ELISA; cross-reacts with human and mouse. Recombinant immunogen mapping to aa 3200-3396. Ideal for extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	Versican Rabbit mAb
Remarks/Alias	WGN; ERVR; GHAP; PG-M; WGN1; CSPG2; Versican
Species	Human
GeneID	1462
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 3200-3396 of human Versican (P13611)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13611
Protein Weight	373kDa
Shipping	Ice bag