

Anti-Versican Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P13611

Item Number: CM0022046

Introduction

Versican rabbit polyclonal antibody validated for WB, IF-P, and ELISA. Detects human, mouse, and rat versican (373 kDa). Ideal for extracellular matrix, cell motility, and differentiation research.

[Learn More](#)

Parameter	Value
Product Name	Versican Rabbit pAb
Remarks/Alias	WGN; ERVR; GHAP; PG-M; WGN1; CSPG2; Versican
Species	Human
GeneID (Human)	1462
GeneID	1462
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1400-1586 of human Versican (NP_004376.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13611
Protein Weight	373kDa
Shipping	Ice bag