

Anti-Vgll3 Rabbit Polyclonal Antibody For Wb And Elisa - A8Mv65

Item Number: CM0016183

Introduction

High-quality rabbit polyclonal antibody targeting human VGLL3 transcription cofactor. Validated for Western blot (WB) and ELISA, with reactivity to rat. Ideal for placental gene regulation studies and assay development.

[Learn More](#)

Parameter	Value
Product Name	VGLL3 Rabbit pAb
Remarks/Alias	VGL3; VGL-3; VGLL3
Species	Human
GenelD (Human)	389136
GenelD	389136
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	A8MV65
Protein Weight	36kDa
Shipping	Ice bag