


Anti-Vgll4 Polyclonal Antibody For Wb, Elisa - Q14135

Item Number: CM0020815



Introduction

High-quality rabbit polyclonal antibody targeting human VGLL4, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying transcriptional regulation and nuclear coactivator mechanisms.

[Learn More](#)

Parameter	Value
Product Name	VGLL4 Rabbit pAb
Remarks/Alias	VGL-4; VGLL4
Species	Human
Gene ID	9686 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 96-296 of human VGLL4 (NP_001121691.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14135
Protein Weight	31kDa
Shipping	Ice bag