


Anti-Vgll4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q80V24

Item Number: CM0019211



Introduction

Anti-Vgll4 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for transcription cofactor Vgll4 research. Recombinant immunogen (aa 1-287).

[Learn More](#)

Parameter	Specification
Product Name	Vgll4 Rabbit pAb
Remarks/Alias	VGL-4, Vgll4
Species	Mouse
GeneID (Human)	9686
GeneID (Mouse)	232334
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-287 of mouse Vgll4 (NP_808351.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q80V24
Protein Weight	31 kDa
Shipping	Ice bag