

Anti-Vglut1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9P2U7

Item Number: CM0020903

Introduction

VGLUT1 Rabbit Polyclonal Antibody for WB, IF-P, ELISA applications. Validated in Human, Mouse, Rat. Targets vesicular glutamate transporter 1 (SLC17A7) for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	VGLUT1 Rabbit pAb
Remarks/Alias	BNPI; VGLUT1; SLC17A7
Species	Human
GeneID (Human)	57030
GeneID	57030
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VGLUT1 (NP_064705.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2U7
Protein Weight	62kDa
Shipping	Ice bag