


Anti-Vhl Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P40337

Item Number: CM0017122



Introduction

Rabbit polyclonal antibody targeting VHL (pVHL, RCA1). Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for HIF pathway and tumor suppression studies.

[Learn More](#)

Parameter	Value
Product Name	VHL Rabbit pAb
Remarks/Alias	RCA1; VHL1; pVHL; HRCA1; VHL
Species	Human
GeneID (Human)	7428
GeneID	7428
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 87-186 of human VHL (NP_000542.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40337
Protein Weight	24kDa
Shipping	Ice bag