


Anti-Biii-Tubulin Monoclonal Antibody For If, Ihc, Elisa - Q13509

Item Number: CM0002273



Introduction

Mouse monoclonal antibody targeting human TUBB3 (β III-Tubulin, Q13509). Verified for IF, IHC-P, and ELISA; cross-reacts with mouse and rat. Essential for studying neuronal development, axon guidance, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	β III-Tubulin Mouse mAb
Remarks/Alias	CDCBM; FEOM3; TUBB4; CDCBM1; CFEOM3; beta-4; CFEOM3A; [KD Validated] β III-Tubulin
Species	Human
GeneID (Human)	10381
GeneID	10381
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 40-260 of human TUBB3 (Q13509).
Source	Mouse
Category	Monoclonal Antibodies
Application	IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13509
Protein Weight	50kDa
Shipping	Ice bag