

Anti-Biii-Tubulin Rabbit Monoclonal Antibody For Wb, If, Ihc-P, Elisa - Q13509

Item Number: CM0003139



Introduction

Rabbit monoclonal antibody targeting human β III-Tubulin (TUBB3), validated for WB, IF, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and oncology research.

[Learn More](#)

Parameter	Value
Product Name	β III-Tubulin Rabbit mAb
Remarks/Alias	CDCBM; FEOM3; TUBB4; CDCBM1; CFEOM3; beta-4; CFEOM3A; β III-Tubulin
Species	Human
GeneID (Human)	10381
GeneID	10381
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13509
Protein Weight	50 kDa
Shipping	Ice bag