

Anti-Biii-Tubulin (Tubb3) Monoclonal Antibody For Wb, Ip, If, Ihc, Elisa - Q13509

Item Number: CM0007312



Introduction

KO Validated anti- β III-Tubulin rabbit monoclonal antibody validated for WB, IP, IF, IHC, ELISA. Cross-reacts with human and mouse. Ideal for axon guidance and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] β III-Tubulin Rabbit mAb
Remarks/Alias	CDCBM, FEOM3, TUBB4, CDCBM1, CFEOM3, beta-4, CFEOM3A, Tubulin beta-3 chain, Tubulin beta-4 chain, Tubulin beta-III
Species	Human
GeneID (Human)	10381
Immunogen	Synthetic peptide corresponding to amino acids 351-450 of β III-Tubulin (NP_006077.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13509
Protein Weight	50kDa
Shipping	Ice bag