

Anti-Biii-Tubulin (Tubb3) Rabbit Polyclonal Antibody For Wb And Elisa - Q13509

Item Number: CM0022127



Introduction

β III-Tubulin Rabbit Polyclonal Antibody validated for WB and ELISA in human, mouse, and rat samples. Ideal for microtubule dynamics, neuronal development, and tumor marker studies. 50kDa - Q13509.

[Learn More](#)

Parameter	Value
Product Name	β III-Tubulin Rabbit pAb
Remarks/Alias	CDCBM; FEOM3; TUBB4; CDCBM1; CFEOM3; beta-4; CFEOM3A; β III-Tubulin
Species	Human
GeneID (human)	10381
GeneID	10381
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 40-260 of human β III-Tubulin/ β 3-Tubulin (NP_006077.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13509
Protein Weight	50kDa
Shipping	Ice bag