

Anti-Villin-1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P09327

Item Number: CM0006757



Introduction

Rabbit monoclonal antibody against human Villin-1, validated for WB, IF-P, IHC-P, and ELISA applications. Reacts with human, mouse, and rat samples. Suitable for studying actin dynamics, epithelial cell biology, and microvilli structure.

[Learn More](#)

Parameter	Value
Product Name	Villin1 Rabbit mAb
Remarks/Alias	VIL; D2S1471; Villin1
Species	Human
GeneID (Human)	7429
GeneID	7429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-827 of human Villin1 (NP_009058.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09327
Protein Weight	93kDa
Shipping	Ice bag