

Anti-Villin-Like Protein (Vill) Polyclonal Antibody For Wb, If/Icc, Elisa - O15195

Item Number: CM0014913

Introduction

CamelBio offers VILL Rabbit pAb, a polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, rat villin-like protein, a potential tumor suppressor. Ideal for cancer research and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	VILL Rabbit pAb
Remarks/Alias	VILL
Species	Human
GeneID (human)	50853
GeneID	50853
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 637-856 of human VILL (NP_056957.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15195
Protein Weight	96kDa
Shipping	Ice bag