

Anti-Villin Monoclonal Antibody For Wb, Ihc-P, Elisa - P09327

Item Number: CM0004774



Introduction

Recombinant rabbit monoclonal antibody targeting human Villin-1 for Western blot, immunohistochemistry, and ELISA. Cross-reacts with mouse and rat, suitable for studying actin cytoskeleton in intestinal and renal epithelium.

[Learn More](#)

Parameter	Value
Product Name	Villin Rabbit mAb
Remarks/Alias	VIL; D2S1471; Villin
Species	Human
Gene ID	7429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-827 of human Villin1 (NP_009058.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09327
Protein Weight	93kDa
Shipping	Ice bag