

Anti-Villin1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P09327

Item Number: CM0003177



Introduction

Monoclonal antibody against Villin1 (VIL; D2S1471) validated for Western blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying actin-binding protein function, microvilli formation, and epithelial cell biology.

[Learn More](#)

Parameter	Value
Product Name	Villin1 Mouse mAb
Remarks/Alias	VIL; D2S1471; Villin1
Species	Human
GeneID	7429
Immunogen	Recombinant protein of human Villin1
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09327
Protein Weight	93kDa
Shipping	Ice bag