


Anti-Vimentin (Phospho S83) Rabbit Polyclonal Antibody For Wb And Elisa - P08670

Item Number: CM0001012



Introduction

Phospho-specific rabbit polyclonal antibody against vimentin at S83, validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for mesenchymal cell, cytoskeleton, and migration research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Vimentin-S83 Rabbit pAb
Remarks/Alias	VIM; CTRCT30; HEL113; vimentin; Phospho-Vimentin-S83
Species	Human
GeneID (Human)	7431
GeneID	7431
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S83 of human VIM (NP_003371.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08670
Protein Weight	54kDa
Shipping	Ice bag