


Anti-Vimentin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08670

Item Number: CM0030107



Introduction

Rabbit polyclonal antibody targeting human, mouse, rat Vimentin, validated for WB, IHC-P, IF/ICC, ELISA. Detects the 54kDa intermediate filament protein involved in cell migration, structural support, and tissue integrity.

[Learn More](#)

Parameter	Value
Product Name	Vimentin Rabbit pAb
Remarks/Alias	CTRCT30; HEL113; Vimentin; VIM; vimentin
Species	Human
Gene ID (Human)	7431
Gene ID	7431
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Vimentin (NP_003371.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08670
Protein Weight	54kDa
Shipping	Ice bag