


Anti-Vimentin Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ip, If/Icc, Ihc-P, Elisa - P08670

Item Number: CM0003625



Introduction

KD-validated rabbit monoclonal antibody against Vimentin, suitable for Western blot, immunoprecipitation, immunofluorescence, and immunohistochemistry. Reacts with human, mouse, and rat Vimentin. Immunogen peptide from 367-466 of human Vimentin.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Vimentin Rabbit mAb
Remarks/Alias	CTRCT30, HEL113, Vimentin, VIM, vimentin, in
Species	Human
Gene ID (Human)	7431
Gene ID	7431
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 367-466 of human Vimentin (P08670)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08670
Protein Weight	54 kDa
Shipping	Ice bag