


Anti-Vinculin [Kd Validated] Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P18206

Item Number: CM0009122



Introduction

Anti-Vinculin rabbit monoclonal antibody for WB, IF-P, IHC-P, and ELISA. Detects human, mouse, rat vinculin. KD validated; suitable for cell adhesion and focal adhesion research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Vinculin Rabbit mAb
Remarks/Alias	MV; MVCL; CMD1W; CMH15; HEL114
Species	Human
GeneID (Human)	7414
GeneID	7414
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-258 of human Vinculin (NP_054706.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18206
Protein Weight	124kDa
Shipping	Ice bag