


Anti-Vinculin Monoclonal Antibody For Wb, If-P, Elisa - P18206

Item Number: CM0011431



Introduction

Anti-Vinculin rabbit monoclonal antibody validated for WB, IF-P, and ELISA. Reacts with human, mouse, and rat Vinculin (Swiss-Prot: P18206). Ideal for studying cell adhesion, focal adhesions, and cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	Vinculin Rabbit mAb
Remarks/Alias	MV; MVCL; CMD1W; CMH15; HEL114; Vinculin
Species	Human
Gene ID (Human)	7414
Gene ID	7414
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 800-900 of human Vinculin (NP_054706.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18206
Protein Weight	124kDa
Shipping	Ice bag