


Anti-Vinculin Monoclonal Antibody For Wb, Ihc-P, Elisa - P18206

Item Number: CM0007505



Introduction

Anti-Vinculin rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat vinculin (P18206). Ideal for cell adhesion, focal adhesion, and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Vinculin Rabbit mAb
Remarks/Alias	MV; MVCL; CMD1W; CMH15; HEL114; Vinculin
Species	Human
GeneID (Human)	7414
GeneID	7414
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18206
Protein Weight	124kDa
Shipping	Ice bag