


Anti-Vinculin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P18206

Item Number: CM0021172



Introduction

High-quality anti-Vinculin rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Specific for human, mouse, rat vinculin (P18206). Ideal for cell adhesion and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Vinculin Rabbit pAb
Remarks/Alias	MV; MVCL; CMD1W; CMH15; HEL114
Species	Human
GeneID (human)	7414
GeneID	7414
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-258 of human Vinculin (NP_054706.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18206
Protein Weight	124kDa
Shipping	Ice bag