


Anti-Vip Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P01282

Item Number: CM0006589



Introduction

High-performance anti-VIP rabbit monoclonal antibody for WB, IF/ICC, IHC-P, ELISA. Detects human VIP, a neuropeptide involved in vasodilation, CNS function, and smooth muscle relaxation. Ideal for neuroscience and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	VIP Rabbit mAb
Remarks/Alias	PHM27; VIP
Species	Human
GeneID (Human)	7432
GeneID	7432
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P,ELISA
Cross Reactivity	Human
SWISS	P01282
Protein Weight	19kDa
Shipping	Ice bag