

Anti-Vipas39 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H9C1

Item Number: CM0027338

Introduction

Rabbit polyclonal antibody against human VIPAS39 (SPE39, VIPAR). Suitable for Western blot (WB) and ELISA. Cross-reacts with human and mouse. Useful for studying endosomal trafficking and epithelial polarity.

[Learn More](#)

Parameter	Value
Product Name	VIPAS39 Rabbit pAb
Remarks/Alias	SPE39; VIPAR; SPE-39; VPS16B; hSPE-39; C14orf133; VIPAS39
Species	Human
GeneID (Human)	63894
GeneID	63894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human VIPAS39 (NP_071350.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H9C1
Protein Weight	57kDa
Shipping	Ice bag