


Anti-Vipr1 Polyclonal Antibody For Wb, Elisa - P32241

Item Number: CM0014721



Introduction

Rabbit polyclonal antibody against VIPR1 (VPAC1 receptor) for Western Blot and ELISA. Targets human VIPR1, a G protein-coupled receptor activated by VIP and PACAP. Immunogen: aa 31-150.

[Learn More](#)

Parameter	Value
Product Name	VIPR1 Rabbit pAb
Remarks/Alias	II; HVR1; RDC1; V1RG; VIPR; VIRG; VAPC1; VPAC1; VPAC1R; VIP-R-1; VPCAP1R; PACAP-R2; PACAP-R-2; VIPR1
Species	Human
GeneID (human)	7433
GeneID	7433
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 31-150 of human VIPR1 (NP_004615.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P32241
Protein Weight	52kDa
Shipping	Ice bag