

Anti-Virma/Kiaa1429 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q69Yn4

Item Number: CM0020760



Introduction

High-quality rabbit polyclonal antibody targeting human VIRMA/KIAA1429, validated for WB, IHC-P, ELISA. Recognizes human, mouse, and rat. Ideal for studying m6A methylation and RNA processing.

[Learn More](#)

Parameter	Value
Product Name	VIRMA/KIAA1429 Rabbit pAb
Remarks/Alias	MSTP054; fSAP121; KIAA1429
Species	Human
Gene ID (Human)	25962
Gene ID	25962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1400-1700 of human KIAA1429 (NP_056311.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q69YN4
Protein Weight	202kDa
Shipping	Ice bag