

Anti-Virma/Kiaa1429 Rabbit Monoclonal Antibody (Kd Validated) For Wb, Ihc-P, Elisa - Q69Yn4

Item Number: CM0009364



Introduction

Rabbit monoclonal antibody against human VIRMA/KIAA1429, KD validated. Suitable for Western blot, IHC-P, and ELISA. Cross-reacts with human and mouse. Target for m6A RNA methylation and RNA processing research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] VIRMA/KIAA1429 Rabbit mAb
Remarks/Alias	MSTP054; fSAP121; KIAA1429
Species	Human
Gene ID (Human)	25962
Gene ID	25962
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 221-320 of human VIRMA/KIAA1429 (NP_056311.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q69YN4
Protein Weight	202kDa
Shipping	Ice bag