


Anti-Vista/B7-H5/GI24 Rabbit Polyclonal Antibody For Wb, IHC-P, Elisa - Q9H7M9

Item Number: CM0020833



Introduction

High-quality rabbit polyclonal antibody targeting human VISTA (B7-H5/GI24/PD-1H), a negative immune checkpoint receptor. Suitable for Western blot, immunohistochemistry (paraffin), and ELISA. Shows cross-reactivity with rat. Supports immuno-oncology and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	VISTA / B7-H5 / GI24 Rabbit pAb
Remarks/Alias	VSIR; B7-H5; B7H5; C10orf54; DD1alpha; GI24; PD-1H; PP2135; SISP1; VISTA; V-set immunoregulatory receptor; VISTA / B7-H5 / GI24
Species	Human
GeneID (human)	64115
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-194 of human VISTA / B7-H5 / GI24 (NP_071436.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H7M9
Protein Weight	33kDa
Shipping	Ice bag