


Anti-Vitamin D-Binding Protein Rabbit Monoclonal Antibody For Wb, Elisa - P02774

Item Number: CM0003669



Introduction

Recombinant rabbit monoclonal antibody targeting human Vitamin D-binding protein (DBP, GC, GcMAF). Suitable for Western blot and ELISA; cross-reacts with mouse and rat. For research on vitamin D transport, actin scavenging, and immune response.

[Learn More](#)

Parameter	Value
Product Name	Vitamin D Binding protein Rabbit mAb
Remarks/Alias	DBP; VDB; GRD3; VDBG; VDBP; GcMAF; DBP/GC; Gc-MAF; DBP-maf; HEL-S-51; Vitamin D Binding protein
Species	Human
GeneID (Human)	2638
GeneID	2638
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 305-474 of human Vitamin D Binding protein (NP_000574.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02774
Protein Weight	53kDa
Shipping	Ice bag