


Anti-Vitamin D Binding Protein Rabbit Polyclonal Antibody For Wb, Elisa - P02774

Item Number: CM0028263



Introduction

High-quality rabbit polyclonal antibody targeting human Vitamin D Binding Protein (DBP/Gc). Validated for WB and ELISA, reacts with human and mouse DBP. Suitable for vitamin D transport and immune function studies.

[Learn More](#)

Parameter	Value
Product Name	Vitamin D Binding protein Rabbit pAb
Remarks/Alias	DBP; VDB; GRD3; VDBG; VDBP; GcMAF; DBP/GC; Gc-MAF; DBP-maf; HEL-S-51; Vitamin D Binding protein
Species	Human
Gene ID (Human)	2638
Gene ID	2638
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 324-493 of human Vitamin D Binding protein (NP_001191236.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02774
Protein Weight	53kDa
Shipping	Ice bag