

Anti-Vitronectin Polyclonal Antibody For Wb, Ihc-P, Elisa - P04004

Item Number: CM0019344



Introduction

Rabbit polyclonal antibody against human Vitronectin, validated for WB, IHC-P, and ELISA; detects human, mouse, and rat proteins. Ideal for cell adhesion and complement pathway research.

[Learn More](#)

Parameter	Specification
Product Name	Vitronectin Rabbit pAb
Remarks/Alias	VN; V75; VNT; Vitronectin
Species	Human
GeneID	7448
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 364-478 of human Vitronectin (NP_000629.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04004
Protein Weight	54 kDa
Shipping	Ice bag