


Anti-Vitronectin (S-Protein, V75) Monoclonal Antibody For Wb, Elisa - P04004

Item Number: CM0011606



Introduction

Rabbit monoclonal antibody targeting human Vitronectin (S-protein, V75). Validated for WB and ELISA. Ideal for investigating cell adhesion, complement inhibition, and extracellular matrix biology.

[Learn More](#)

Parameter	Value
Product Name	Vitronectin Rabbit mAb
Remarks/Alias	VN; V75; VNT; Vitronectin
Species	Human
GeneID	7448
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 379-478 of human Vitronectin (P04004).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04004
Protein Weight	54kDa
Shipping	Ice bag