

Anti-Vnn1 Rabbit Polyclonal Antibody For Wb And Elisa - 095497

Item Number: CM0028438



Introduction

Rabbit polyclonal antibody targeting human VNN1 (Pantetheinase), suitable for WB and ELISA applications, with cross-reactivity to mouse. Useful for studying pantothenic acid recycling and inflammatory pathways.

[Learn More](#)

Parameter	Value
Product Name	VNN1 Rabbit pAb
Remarks/Alias	HDLCQ8; Tiff66; VNN1
Species	Human
GeneID (human)	8876
GeneID	8876
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-415 of human VNN1 (NP_004657.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	095497
Protein Weight	57kDa
Shipping	Ice bag