


Anti-Von Willebrand Factor (Vwf) Rabbit Monoclonal Antibody For Ihc-P, Elisa - P04275

Item Number: CM0007546



Introduction

Rabbit monoclonal antibody targeting human von Willebrand factor (VWF) for IHC-P and ELISA applications. Cross-reacts with human and mouse; immunogen spans amino acids 2645-2813. Supports hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	von Willebrand factor (VWF) Rabbit mAb
Remarks/Alias	VWD; F8VWF; von Willebrand factor (VWF)
Species	Human
GeneID (Human)	7450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2645-2813 of human von Willebrand factor (VWF) (NP_000543.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04275
Protein Weight	309kDa
Shipping	Ice bag