


Anti-Von Willebrand Factor (Vwf) Rabbit Monoclonal Antibody For Wb - P04275

Item Number: CM0003056



Introduction

This rabbit monoclonal antibody against human von Willebrand factor (VWF) is validated for WB and ELISA, cross-reacts with mouse. Designed for hemostasis research, it targets the VWF protein involved in platelet adhesion and coagulation factor VIII stabilization.

[Learn More](#)

Parameter	Value
Product Name	von Willebrand factor (VWF) Rabbit mAb
Remarks/Alias	VWD; F8VWF; von Willebrand factor (VWF)
Species	Human
GeneID (Human)	7450
GeneID	7450
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P04275
Protein Weight	309kDa
Shipping	Ice bag